

ISOLATION – BI-DIRECTIONAL

Limited personnel access. Use immune deficient mice in isolation first.

Prior to room entry:



Inside room:



Facility entrance:

1. Shoe covers
2. Gown

Inside room:

1. Sleeves or apron
2. Second pair of gloves
3. Spray gloves & sleeves with Oxivir & wait 5 minutes

Outside room:

1. Bouffant
2. Surgical mask
3. Push thumbs through sleeves of gown
4. Gloves

Depart in your PPE. If additional consumable supplies, equipment, reagents or biologics are needed during study, they must be placed in nonporous secondary containers and either sterilized by autoclaving or vaporized hydrogen peroxide, or disinfected by saturate spraying with Oxivir TB. Containers are then delivered using a **two person** transfer procedure. If mice depart and return to isolation, a **two person** transfer procedure must be used. If reentry is required, don a **Tyvek suit** and follow entry procedures above.

